



PDE11A Antibody, Biotin conjugated

Product Code	CSB-PA875726LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCR9
Immunogen	Recombinant Human Dual 3',5'-cyclic-AMP and -GMP phosphodiesterase 11A protein (42-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual 3',5'-cyclic-AMP and -GMP phosphodiesterase 11A (EC 3.1.4.35) (EC 3.1.4.53) (cAMP and cGMP phosphodiesterase 11A), PDE11A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE11A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.