



RBM38 Antibody, Biotin conjugated

Product Code	CSB-PA880939LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0Z9
Immunogen	Recombinant Human RNA-binding protein 38 protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 38 (CLL-associated antigen KW-5) (HSRNASEB) (RNA-binding motif protein 38) (RNA-binding region-containing protein 1) (ssDNA-binding protein SEB4), RBM38, RNPC1 SEB4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RBM38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.