



HAO2 Antibody, Biotin conjugated

Product Code	CSB-PA882145LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYQ3
Immunogen	Recombinant Human Hydroxyacid oxidase 2 protein (2-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxyacid oxidase 2 (HAOX2) (EC 1.1.3.15) ((S)-2-hydroxy-acid oxidase, peroxisomal) (Cell growth-inhibiting gene 16 protein) (Long chain alpha-hydroxy acid oxidase) (Long-chain L-2-hydroxy acid oxidase), HAO2, HAOX2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HAO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.