




GGA3 Antibody

Product Code	CSB-PA882151LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ52
Immunogen	Recombinant Human ADP-ribosylation factor-binding protein GGA3 protein (342-503AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor-binding protein GGA3 (Golgi-localized, gamma ear-containing, ARF-binding protein 3), GGA3, KIAA0154
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GGA3

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA882151LA01HU at dilution of 1:100