



OLFML3 Antibody, Biotin conjugated

Product Code	CSB-PA885708LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRN5
Immunogen	Recombinant Human Olfactomedin-like protein 3 protein (201-406AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactomedin-like protein 3 (HNOEL-iso) (hOLF44), OLFML3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OLFML3