



ACAD9 Antibody, FITC conjugated

Product Code	CSB-PA887999EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H845
Immunogen	Recombinant Human Acyl-CoA dehydrogenase family member 9, mitochondrial protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-CoA dehydrogenase family member 9, mitochondrial (ACAD-9) (EC 1.3.99.-), ACAD9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACAD9