



SLC6A13 Antibody, FITC conjugated

Product Code	CSB-PA889104LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSD5
Immunogen	Recombinant Human Sodium- and chloride-dependent GABA transporter 2 protein (1-40AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium- and chloride-dependent GABA transporter 2 (GAT-2) (Solute carrier family 6 member 13), SLC6A13, GAT2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC6A13