



VPS29 Antibody, FITC conjugated

Product Code	CSB-PA890661LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBQ0
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 29 protein (68-145AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 29 (hVPS29) (PEP11 homolog) (Vesicle protein sorting 29), VPS29
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS29