



VPS29 Antibody, Biotin conjugated

Product Code	CSB-PA890661LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBQ0
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 29 protein (68-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 29 (hVPS29) (PEP11 homolog) (Vesicle protein sorting 29), VPS29
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.