




B4GALT7 Antibody

Product Code	CSB-PA890666LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBV7
Immunogen	Recombinant Human Beta-1,4-galactosyltransferase 7 protein (52-327AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-1,4-galactosyltransferase 7 (Beta-1,4-GalTase 7) (Beta4Gal-T7) (b4Gal-T7) (EC 2.4.1.-) (UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 7) (UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 7) [Includes: Xylosylprotein 4-beta-galactosyltransferase (EC 2.4.1.133) (Proteoglycan UDP-galactose:beta-xylose beta1,4-galactosyltransferase I) (UDP-galactose:beta-xylose beta-1,4-galactosyltransferase) (XGPT) (XGalT-1) (Xylosylprotein beta-1,4-galactosyltransferase)], B4GALT7, XGALT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	B4GALT7

Image



Western Blot

Positive WB detected in: Mouse brain tissue

All lanes: B4GALT7 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38 kDa

Observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.