




CDC42EP3 Antibody

Product Code	CSB-PA890918LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKI2
Immunogen	Recombinant Human Cdc42 effector protein 3 protein (1-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cdc42 effector protein 3 (Binder of Rho GTPases 2) (MSE55-related Cdc42-binding protein), CDC42EP3, BORG2 CEP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CDC42EP3

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA890918LA01HU at dilution of 1:100