



SHOC2 Antibody, FITC conjugated

Product Code	CSB-PA891468LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ13
Immunogen	Recombinant Human Leucine-rich repeat protein SHOC-2 protein (3-89AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat protein SHOC-2 (Protein soc-2 homolog) (Protein sur-8 homolog), SHOC2, KIAA0862
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SHOC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.