



USP18 Antibody, Biotin conjugated

Product Code	CSB-PA891559LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMW8
Immunogen	Recombinant Human Ubl carboxyl-terminal hydrolase 18 protein (18-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubl carboxyl-terminal hydrolase 18 (EC 3.4.19.-) (43 kDa ISG15-specific protease) (hUBP43) (ISG15-specific-processing protease) (Ubl thioesterase 18), USP18, ISG43
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.