



FZD1 Antibody, Biotin conjugated

Product Code	CSB-PA891570LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UP38
Immunogen	Recombinant Human Frizzled-1 protein (73-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Frizzled-1 (Fz-1) (hFz1) (FzE1), FZD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FZD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.