



RIMS2 Antibody, HRP conjugated

Product Code	CSB-PA892173LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ26
Immunogen	Recombinant Human Regulating synaptic membrane exocytosis protein 2 protein (667-943AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulating synaptic membrane exocytosis protein 2 (Rab-3-interacting molecule 2) (RIM 2) (Rab-3-interacting protein 3), RIMS2, KIAA0751 RAB3IP3 RIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RIMS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.