



RIMS2 Antibody, FITC conjugated

Product Code	CSB-PA892173LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ26
Immunogen	Recombinant Human Regulating synaptic membrane exocytosis protein 2 protein (667-943AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulating synaptic membrane exocytosis protein 2 (Rab-3-interacting molecule 2) (RIM 2) (Rab-3-interacting protein 3), RIMS2, KIAA0751 RAB3IP3 RIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RIMS2