



GAB2 Antibody, FITC conjugated

Product Code	CSB-PA892175LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQC2
Immunogen	Recombinant Human GRB2-associated-binding protein 2 protein (273-484AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GRB2-associated-binding protein 2 (GRB2-associated binder 2) (Growth factor receptor bound protein 2-associated protein 2) (pp100), GAB2, KIAA0571
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GAB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.