





SLC12A7 Antibody

Product Code	CSB-PA893907LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y666
Immunogen	Recombinant Human Solute carrier family 12 member 7 protein (11-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 7 (Electroneutral potassium-chloride cotransporter 4) (K-Cl cotransporter 4), SLC12A7, KCC4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC12A7


Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA893907LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA893907LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA893907LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.