



SLC12A7 Antibody, Biotin conjugated

Product Code	CSB-PA893907LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y666
Immunogen	Recombinant Human Solute carrier family 12 member 7 protein (11-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 7 (Electroneutral potassium-chloride cotransporter 4) (K-Cl cotransporter 4), SLC12A7, KCC4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC12A7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.