



mib1 Antibody, Biotin conjugated

Product Code	CSB-PA896152LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	H9GYU5
Immunogen	Recombinant Zebrafish E3 ubiquitin-protein ligase mib1 protein (1-274AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	mib1
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	mib1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.