



HIF3A Antibody, Biotin conjugated

Product Code	CSB-PA896696LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2N7
Immunogen	Recombinant Human Hypoxia-inducible factor 3-alpha protein (516-669AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hypoxia-inducible factor 3-alpha (HIF-3-alpha) (HIF3-alpha) (Basic-helix-loop-helix-PAS protein MOP7) (Class E basic helix-loop-helix protein 17) (bHLHe17) (HIF3-alpha-1) (Inhibitory PAS domain protein) (IPAS) (Member of PAS protein 7) (PAS domain-containing protein 7), HIF3A, BHLHE17 MOP7 PASD7
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HIF3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.