




# KLF12 Antibody

<b>Product Code</b>	CSB-PA897495LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y4X4
<b>Immunogen</b>	Recombinant Human Krueppel-like factor 12 protein (36-313AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Krueppel-like factor 12 (Transcriptional repressor AP-2rep), KLF12, AP2REP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	KLF12

## Image



### Western Blot

Positive WB detected in: Mouse kidney tissue, Rat liver tissue, Rat brain tissue

All lanes: KLF12 antibody at 2.5µg/ml

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45, 29, 34 kDa

Observed band size: 45 kDa