



KLF12 Antibody, FITC conjugated

Product Code	CSB-PA897495LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4X4
Immunogen	Recombinant Human Krueppel-like factor 12 protein (36-313AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Krueppel-like factor 12 (Transcriptional repressor AP-2rep), KLF12, AP2REP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KLF12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.