



ETV7 Antibody, FITC conjugated

Product Code	CSB-PA897576LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y603
Immunogen	Recombinant Human Transcription factor ETV7 protein (134-218AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor ETV7 (ETS translocation variant 7) (ETS-related protein Tel2) (Tel-related Ets factor) (Transcription factor Tel-2), ETV7, TEL2 TELB TREF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ETV7