



# Goat Anti-Rabbit IgG(H+L) Antibody;HRP-conjugated

<b>Product Code</b>	CSB-PA001501GA01RB
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Goat
<b>Species Reactivity</b>	Rabbit
<b>Tested Applications</b>	ELISA,WB;Recommended dilution:WB:1:2000-1:20000
<b>Form</b>	Liquid
<b>Storage Buffer</b>	0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.
<b>Purification Method</b>	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Product Type</b>	Secondary Antibody
<b>Immunogen Species</b>	Rabbit
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.