



Goat Anti-Rabbit IgG(H+L) Antibody;HRP-conjugated

Product Code	CSB-PA001501GA01RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Goat
Species Reactivity	Rabbit
Tested Applications	ELISA,WB;Recommended dilution:WB:1:2000-1:20000
Form	Liquid
Storage Buffer	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
Purification Method	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Product Type	Secondary Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)