





DUX4 Antibody

Product Code	CSB-PA890667LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX2
Immunogen	Recombinant Human Double homeobox protein 4 protein (217-402AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Double homeobox protein 4 (Double homeobox protein 10), DUX4, DUX10
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	DUX4

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA890667LA01HU at dilution of 1: 100



Western Blot

Positive WB detected in: Hela whole cell lysate, PC-3 whole cell lysate, A549 whole cell lysate, HepG2 whole cell lysate

All lanes: DUX4 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.