




PRSS1 Antibody

Product Code	CSB-PA12309A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07477
Immunogen	Recombinant Human Trypsin-1 protein (104-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trypsin-1 (EC 3.4.21.4) (Beta-trypsin) (Cationic trypsinogen) (Serine protease 1) (Trypsin I) [Cleaved into: Alpha-trypsin chain 1; Alpha-trypsin chain 2], PRSS1, TRP1 TRY1 TRYP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRSS1

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA12309A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.