




DSG3 Antibody

Product Code	CSB-PA007205YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32926
Immunogen	Recombinant Human Desmoglein-3 protein (50-615AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmoglein-3 (130 kDa pemphigus vulgaris antigen) (PVA) (Cadherin family member 6), DSG3, CDHF6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSG3

Image



Western Blot

Positive WB detected in: Mouse brain tissue

All lanes: DSG3 antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 108 kDa

observed band size: 108 kDa