




Ripk1 Antibody

Product Code	CSB-PA720181LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q60855
Immunogen	Recombinant Mouse Receptor-interacting serine/threonine-protein kinase 1 protein
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Ripk1

Image



Western Blot

Positive WB detected in: Mouse Kidney tissue, HeLa whole cell lysate

All lanes: Ripk1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75 kDa

Observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.