



DUX4 Antibody, HRP conjugated

Product Code	CSB-PA890667LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX2
Immunogen	Recombinant Human Double homeobox protein 4 protein (217-402AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Double homeobox protein 4 (Double homeobox protein 10), DUX4, DUX10
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	DUX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.