



DUX4 Antibody, Biotin conjugated

Product Code	CSB-PA890667LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX2
Immunogen	Recombinant Human Double homeobox protein 4 protein (217-402AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Double homeobox protein 4 (Double homeobox protein 10), DUX4, DUX10
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	DUX4