



DSG3 Antibody, Biotin conjugated

Product Code	CSB-PA007205YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32926
Immunogen	Recombinant Human Desmoglein-3 protein (50-615AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmoglein-3 (130 kDa pemphigus vulgaris antigen) (PVA) (Cadherin family member 6), DSG3, CDHF6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSG3