



# Ripk1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA720181LC01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q60855
<b>Immunogen</b>	Recombinant Mouse Receptor-interacting serine/threonine-protein kinase 1 protein (1-656AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP) (Receptor-interacting protein 1) (RIP-1) (Serine/threonine-protein kinase RIP), Ripk1, Rinp Rip
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	Ripk1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.