



Ripk1 Antibody, Biotin conjugated

Product Code	CSB-PA720181LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q60855
Immunogen	Recombinant Mouse Receptor-interacting serine/threonine-protein kinase 1 protein (1-656AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-interacting serine/threonine-protein kinase 1 (EC 2.7.11.1) (Cell death protein RIP) (Receptor-interacting protein 1) (RIP-1) (Serine/threonine-protein kinase RIP), Ripk1, Rinp Rip
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Ripk1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.