



CCR3 Antibody, HRP conjugated

Product Code	CSB-PA004842LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51677
Immunogen	Recombinant Human C-C chemokine receptor type 3 protein (306-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C chemokine receptor type 3 (C-C CKR-3) (CC-CKR-3) (CCR-3) (CCR3) (CKR3) (Eosinophil eotaxin receptor) (CD antigen CD193), CCR3, CMKBR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.