



CCR3 Antibody, HRP conjugated

Product Code	CSB-PA004842LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51677
Immunogen	Recombinant Human C-C chemokine receptor type 3 protein (306-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C chemokine receptor type 3 (C-C CKR-3) (CC-CKR-3) (CCR-3) (CCR3) (CKR3) (Eosinophil eotaxin receptor) (CD antigen CD193), CCR3, CMKBR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCR3