



CBLL2 Antibody, FITC conjugated

Product Code	CSB-PA836684LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N7E2
Immunogen	Recombinant Human E3 ubiquitin-protein ligase CBLL2 protein (200-413AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase ZNF645 (EC 2.3.2.27) (RING-type E3 ubiquitin transferase ZNF645) (Zinc finger protein 645), ZNF645
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF645
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.