



TRAF3 Antibody, Biotin conjugated

Product Code	CSB-PA024148LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13114
Immunogen	Recombinant Human TNF receptor-associated factor 3 protein (1-40AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TNF receptor-associated factor 3 (EC 2.3.2.27) (CAP-1) (CD40 receptor-associated factor 1) (CRAF1) (CD40-binding protein) (CD40BP) (LMP1-associated protein 1) (LAP1) (RING-type E3 ubiquitin transferase TRAF3), TRAF3, CAP1 CRAF1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TRAF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.