



Recombinant Mouse 5-hydroxytryptamine receptor 2A (Htr2a), partial

Product Code	CSB-BP010887MO1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35363
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	KTYRSAFSRYIQCQYKENRKPLQLILVNTIPTLAYKSSQLQVGQKKNQEDAEP TANDCS MVTLGNQHSEEMCTDNIETVNEKVSCV
Source	Baculovirus
Target Names	Htr2a
Protein Names	Recommended name: 5-hydroxytryptamine receptor 2A Short name= 5-HT-2 Short name= 5-HT-2AAlternative name(s): Serotonin receptor 2A
Expression Region	385-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial, the last Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.