



Recombinant Human UDP-glucuronosyltransferase 2B7 (UGT2B7), partial

Product Code	CSB-EP025598HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16662
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GKVLVWAAEYSHWMNIKTILDELIQRGHEVTVLASSASILFDPNNSALKIEIYP TSLTK TELENFIMQQIKRWSDLPKDTFWLYFSQVQEIMSFIDITRKFCCKDVVSNKKFM KKVQES RFDVIFADAIFFPCSELLAELFNIPFVYSLSFSPGYTFEKHSGGFIFPPSYVPVVM SELT QMTFMERVKNMIYVLYFDFWFEIFDMKKWDQFYSEVLGRPTTLSETMGKADV WLIRNSWN FQFPHLLPNVDFVGGHLHCKPAKPLPKEMEDFVQSSGENGVVVFSLGSMVSN MTEERANV IASALAQIPQKVLWRFDGNKPDTLGLNTRLYKWIPQNDLLGHPKTRAFITHGGA NGIYEA IYHGIPMVGIPLFADQPDNIAHMKARGAAVRVDFNTMSSTDLLNALKRVINDPS YKENVM KLSRIQHDQPVKPLDRAVFWIEFVMRHKGAKHLRVAHDLTWFQYHSLD
Source	E.coli
Target Names	UGT2B7
Protein Names	Recommended name: UDP-glucuronosyltransferase 2B7 Short name= UDPGT 2B7 EC= 2.4.1.17 Alternative name(s): 3,4-catechol estrogen-specific UDPGT UDP-glucuronosyltransferase 2B9 Short name= UDPGT 2B9 UDPGTh-2
Expression Region	24-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.