



# Recombinant Human UDP-glucuronosyltransferase 2B7 (UGT2B7), partial

<b>Product Code</b>	CSB-YP025598HU1
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P16662
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GKVLVWAAEYSHWMNIKTILDELIQRGHEVTVLASSASILFDPNNSALKIEIYP TSLTK TELENFIMQQIKRWSDLPKDTFWLYFSQVQEIMSFIDITRKFCCKDVVSNKKFM KKVQES RFDVIFADAIFFPCSELLAELFNIPFVYSLSFSPGYTFEKHSGGFIFPPSYVPVVM SELT QMTFMERVKNMIYVLYFDFWFEIFDMKKWDQFYSEVLGRPTTLSETMGKADV WLIRNSWN FQFPHLLPNVDFVGGHLHCKPAKPLPKEMEDFVQSSGENGVVVFSLGSMVSN MTEERANV IASALAQIPQKVLWRFDGNKPDTLGLNTRLYKWIPQNDLLGHPKTRAFITHGGA NGIYEA IYHGIPMVGIPLFADQPDNIAHMKARGAAVRVDFNTMSSTDLLNALKRVINDPS YKENVM KLSRIQHDQPVKPLDRAVFWIEFVMRHKGAKHLRVAHDLTWFQYHSLD
<b>Source</b>	Yeast
<b>Target Names</b>	UGT2B7
<b>Protein Names</b>	Recommended name: UDP-glucuronosyltransferase 2B7 Short name= UDPGT 2B7 EC= 2.4.1.17 Alternative name(s): 3,4-catechol estrogen-specific UDPGT UDP-glucuronosyltransferase 2B9 Short name= UDPGT 2B9 UDPGTh-2
<b>Expression Region</b>	24-429
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.