



DPM1 Antibody, HRP conjugated

Product Code	CSB-PA007134LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60762
Immunogen	Recombinant Human Dolichol-phosphate mannosyltransferase subunit 1 protein (2-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dolichol-phosphate mannosyltransferase subunit 1 (EC 2.4.1.83) (Dolichol-phosphate mannose synthase subunit 1) (DPM synthase subunit 1) (Dolichyl-phosphate beta-D-mannosyltransferase subunit 1) (Mannose-P-dolichol synthase subunit 1) (MPD synthase subunit 1), DPM1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	DPM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.