




ETV1 Antibody

Product Code	CSB-PA007854LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50549
Immunogen	Recombinant Human ETS translocation variant 1 protein (61-332AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ETS translocation variant 1 (Ets-related protein 81), ETV1, ER81
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ETV1

Image



Western Blot

Positive WB detected in: Mouse liver tissue, Rat liver tissue

All lanes: ETV1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56, 54, 53, 44, 49, 51 kDa

Observed band size: 56 kDa