




# ETV1 Antibody

<b>Product Code</b>	CSB-PA007854LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P50549
<b>Immunogen</b>	Recombinant Human ETS translocation variant 1 protein (61-332AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	ETS translocation variant 1 (Ets-related protein 81), ETV1, ER81
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	ETV1

## Image



### Western Blot

Positive WB detected in: Mouse liver tissue, Rat liver tissue

All lanes: ETV1 antibody at 3µg/ml

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56, 54, 53, 44, 49, 51 kDa

Observed band size: 56 kDa