



FPR3 Antibody, FITC conjugated

Product Code	CSB-PA008856LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25089
Immunogen	Recombinant Human N-formyl peptide receptor 3 protein (307-353AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-formyl peptide receptor 3 (FMLP-related receptor II) (FMLP-R-II) (Formyl peptide receptor-like 2), FPR3, FPRH1 FPRL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FPR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.