



IFITM3 Antibody, Biotin conjugated

Product Code	CSB-PA011026LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01628
Immunogen	Recombinant Human Interferon-induced transmembrane protein 3 protein (1-57AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced transmembrane protein 3 (Dispanin subfamily A member 2b) (DSPA2b) (Interferon-inducible protein 1-8U), IFITM3
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	IFITM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.