



KIR2DS2 Antibody, Biotin conjugated

Product Code	CSB-PA012360LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43631
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DS2 protein (104-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell immunoglobulin-like receptor 2DS2 (CD158 antigen-like family member J) (MHC class I NK cell receptor) (NK receptor 183 Act1) (Natural killer-associated transcript 5) (NKAT-5) (p58 natural killer cell receptor clone CL-49) (p58 NK receptor CL-49) (CD antigen CD158j), KIR2DS2, CD158J NKAT5
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KIR2DS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.