



PIGA Antibody, Biotin conjugated

Product Code	CSB-PA017965LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37287
Immunogen	Recombinant Human Phosphatidylinositol N-acetylglucosaminyltransferase subunit A protein (443-484AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol N-acetylglucosaminyltransferase subunit A (EC 2.4.1.198) (GlcNAc-PI synthesis protein) (Phosphatidylinositol-glycan biosynthesis class A protein) (PIG-A), PIGA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PIGA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.