



RGL2 Antibody, HRP conjugated

Product Code	CSB-PA019625LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15211
Immunogen	Recombinant Human Ral guanine nucleotide dissociation stimulator-like 2 protein (145-328AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ral guanine nucleotide dissociation stimulator-like 2 (RalGDS-like 2) (RalGDS-like factor) (Ras-associated protein RAB2L), RGL2, RAB2L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.